

# CHOICE FOR FREEDOM



## WIDE RANGE OF DISPOSABLES FOR 12H OF USE


- Use a new secufill®/secufill® XL patient line for each patient
- Use with prefilled bag, bag to fill or both
- Use with vials or bags
- Mono or bi-liquid

### manyFlow® S softbag



Transfer line for soft bag filling with CM\* in ScanBag®

### manyFlow® S ScanBag®



Injection line for CM\* in ScanBag®

### manyFlow® S vial



Transfer line for CM\* in vial



### manyFlow® D CM softbag / NaCl vial



Transfer line for soft bag filling with CM\* in ScanBag® and NaCl in vial

### manyFlow® D CM ScanBag® / NaCl vial



Transfer line for CM\* in ScanBag® and NaCl in vial

### manyFlow® D 2 vial



Transfer line for CM\* in vial and NaCl in vial

### secufill®



Single use safety patient line with double-level valve - 23cm

### secufill® XL



Single use safety patient line with double-level valve - 100cm

\*CM = contrast media